

## Outcomes of Mechanically Ventilated COVID-19 Patients



Millions of people worldwide have been infected with COVID-19, with many requiring hospitalisation. Approximately 30% of hospitalised patients with COVID-19 are started on invasive ventilation for severe pneumonia. Invasive ventilation portends a poor prognosis as approximately 40% of these patients do not survive hospitalisation. Even those who survive usually have a prolonged duration of ventilatory support. When the clinical condition of such patients is stabilized, they are usually transferred to long-term acute care hospitals.

A study was conducted to describe the clinical characteristics and outcomes of adult patients with COVID-19 requiring weaning from prolonged mechanical ventilation. The study included patients admitted to two long-term acute care hospitals in the U.S. from April 1, 2020 to March 31, 2020.

During the study period, 215 COVID-19 patients were transferred to long-term acute care hospitals. Of these, 158 patients had respiratory failure caused by SARS-CoV-2 pneumonia and required weaning from prolonged mechanical ventilation. Prior to long-term acute care transfer, the median length of stay at the acute care hospital was 41 days and the median number of ventilator days was 35. The median age of patients included in the study was 69 years. 34.8% of the patients were women and 91.8% had at least one comorbidity (primarily hypertension and diabetes).

The percent of weaning success was 70.9%. Success weaning duration was 8 days and mortality was 9.6%. 19% of patients were discharged home; 70.3% were discharged to other facilities while 1.3% remained in the long-term acute care hospitals.

Overall, this analysis shows that patients with COVID-19 who were transferred to long-term acute care were successfully weaned from prolonged mechanical ventilation.

Source: [Critical Care Medicine](#)

Image Credit: iStock

Published on : Tue, 17 Aug 2021