



---

## Schiller Launches CARDIOVIT AT-1 G2 at MEDICA 2016



---

### The power of simplicity CARDIOVIT AT-1 G2

SCHILLER is very pleased to introduce the 2nd generation of its most successful ECG device: Building on the solidly proven AT-1 and enhanced with the latest technology, the AT-1 G2 distinguishes itself by its outstanding signal quality and the newest interpretation algorithm. User friendliness is guaranteed with step-by-step workflow and easy patient data entry. A colour screen and an easy-to-clean keyboard complete this intuitive and reliable electrocardiograph.

#### State of the Art

Thanks to a very high sampling frequency and a large frequency band, the CARDIOVIT AT-1 G2 provides optimal signal quality for adult and paediatric ECG. The high-resolution thermal printer is compatible with roll and z-folded paper, accommodating customer preference. Connectivity via USB is also available for software updates.

□

#### Ergonomic Design

The CARDIOVIT AT-1 G2 has been designed with ease-of-use in mind: patient data can be entered at the touch of a button; with a 5" colour display for ECG preview and simplified keyboard, taking an ECG has never been simpler.

#### Step by Step Workflow

1. Patient data entry (at the touch of a button)
2. Real-time signal preview and quality check, with lead reversal detection
3. Post-acquisition review

### **Advanced Clinical Interpretation**

In addition to signal quality check and lead reversal detection, the renowned 12-lead ECG interpretation program ETM is available.

For more information, come and join us at the MEDICA in Dusseldorf (hall 9, stand E05).

We look forward to your visit!

**Source & Image Credit:** [Schiller](#)

Published on : Mon, 14 Nov 2016